

PANDAA[®] EBOV Qualification Panel

Instructions for Use

PANDAA[®] EBOV Qualification Panel is for Research Use Only (RUO)

NOT FOR RESALE

RUO

For research use only

REF

4111096



-15 °C to -25 °C



Contents sufficient for 10 × 20 reactions



Aldatu Biosciences, Inc.
313 Pleasant Street
Watertown, MA 02472
USA

Overview	4
Intended Use (RUO)	4
Contents and Storage	4
Product Description	4
Storage and Stability	5
PANDAA EBOV Qualification Panel Step-by-Step Instructions	5
Before Beginning	5
Running the Qualification Panel	5
Real-Time PCR Protocol	7
Data Analysis	7
Customer and Technical Support	8
Disclaimers	8
Explanation of Symbols	9

Overview

Intended Use (RUO)

The PANDAA EBOV Qualification Panel is a panel of synthetic, non-infectious RNA for use exclusively with the PANDAA® Ebola (Aldatu Biosciences, catalog # 4011096) *in vitro* real-time RT-PCR assays for the amplification and pan-species detection of ebolavirus RNA. This Qualification Panel may be used for training technical personnel in the PANDAA® Ebola assay, internal proficiency testing as a component of a quality assurance program, and for qualifying real-time PCR instruments for use with the PANDAA® Ebola assay.

For Research Use Only. Not for use in diagnostic procedures.

Contents and Storage

Product Description

The PANDAA EBOV Qualification Panel includes ten (10) samples **made from synthetic, non-infectious RNA**. Nine (9) positive panel members covering the target regions for PANDAA assays in the genome of seven ebolavirus isolates. The negative panel member (sample G) consists of the synthetic, non-infectious RNA from the equivalent genomic region of Marburg virus a closely related filovirus. Synthetic RNA is prepared in a background of Internal Control RNA and human genomic DNA.

Each Qualification Panel sample is provided at 200 µL and 10 µL is used in a 20 µL PANDAA reaction.

Label	Virus	Species	Isolate	Conc.	Copies/Reaction
A	Bundibugyo virus	<i>Orthoebolavirus bundibugyoense</i>	Butalya-811250	5 cp/µL	50 cp/rxn
B	Reston virus	<i>Orthoebolavirus restonense</i>	Pennsylvania	5 cp/µL	50 cp/rxn
C	Sudan virus	<i>Orthoebolavirus sudanense</i>	Boniface	5 cp/µL	50 cp/rxn
D	Tai Forest virus	<i>Orthoebolavirus taiense</i>	Pauléoula-CI	5 cp/µL	50 cp/rxn
E	Ebola virus	<i>Orthoebolavirus zairense</i>	Makona	5 cp/µL	50 cp/rxn
F	Sudan virus	<i>Orthoebolavirus sudanense</i>	Gulu	5 cp/µL	50 cp/rxn
G	Marburg virus *	<i>Orthomarburgvirus marburgense</i>	Musoke	500 cp/µL	5,000 cp/rxn
H	Ebola virus	<i>Orthoebolavirus zairense</i>	Mayinga	5,000 cp/µL	50,000 cp/rxn
I	Ebola virus	<i>Orthoebolavirus zairense</i>	Mayinga	500 cp/µL	5,000 cp/rxn
J	Ebola virus	<i>Orthoebolavirus zairense</i>	Mayinga	5 cp/µL	50 cp/rxn

* **Note:** Marburg virus is a filovirus closely related to ebolaviruses and acts as a specificity negative control for the PANDAA EBOV Qualification Panel.

Storage and Stability

The PANDAA® EBOV Qualification Panel is shipped on dry ice, and the kit components will arrive frozen. All components should be stored between -15 °C and -25 °C upon arrival. After thawing, the Qualification Panel samples can be stored between 2 °C and 8 °C for up to 72 hours. More than two freeze-thaw cycles should be avoided to ensure preserve sample integrity.

PANDAA EBOV Qualification Panel Step-by-Step Instructions

Before Beginning

Follow all steps in the instructions for use for the PANDAA® Ebola assay to prepare the PANDAA mastermix.

Running the Qualification Panel

Run each Qualification Panel sample in four replicates alongside four replicates each of the PANDAA® Ebola assay kit Positive and Negative Controls.

In total, you will run 48 samples (12 × 4), which uses one half of a PANDAA assay kit.

- 1. Prepare samples:** Spin the sample tubes briefly to collect drops that may have formed on the sides of the tubes. Gently vortex the sample tubes and spin briefly again.
- 2. Add sample or controls:** Add 10 µL Qualification Panel sample, Positive Control, or Negative Control into the top of the mastermix in each well / tube. **Do not pipette up and down to mix.**
Refer to the sample loading plate maps on page 6.
- 3. Seal and spin (plate only):** Seal the plate with optical adhesive film or cap the reaction tubes. If using a 96-well plate, spin the plate briefly before loading into the real-time PCR instrument.

Qualification Panel Sample Loading Plate Map (96-well Plate)

Load samples in the configuration shown in Figure 1.

	1	2	3	4	5	6	7	8	9	10	11	12
A	A	B	C	D	E	F	G	H	I	J	POS	NEG
B	A	B	C	D	E	F	G	H	I	J	POS	NEG
C	A	B	C	D	E	F	G	H	I	J	POS	NEG
D	A	B	C	D	E	F	G	H	I	J	POS	NEG
E												
F												
G												
H												

Figure 1: Loading map for 96-well plate.

Qualification Panel Sample Loading Map (Optical Reaction Tubes)

Mic optical reaction tubes are pre-packed into a loading rack. Follow the loading rack map in Figure 2 for mastermix and sample loading.

Follow Table 1 for sample loading for the Mic (Bio Molecular Systems) and Rotor-Gene Q (Qiagen). Sample A should be loaded into rotor positions 1–4 with samples B through J then loaded sequentially. The Positive Control should be loaded in positions 41–44 and the Negative Control in positions 45–48.

	1	2	3	4	5	6
A	A	C	E	G	I	POS
B	A	C	E	G	I	POS
C	A	C	E	G	I	POS
D	A	C	E	G	I	POS
E	B	D	F	H	J	NEG
F	B	D	F	H	J	NEG
G	B	D	F	H	J	NEG
H	B	D	F	H	J	NEG

Figure 2: Loading rack map for Mic reaction tubes.

Sample	Rotor Position
A	1–4
B	5–8
C	9–12
D	13–16
E	17–20
F	21–24
G	25–28
H	29–32
I	33–36
J	37–40
POS	41–44
NEG	45–48

Table 1: Rotor loading

Real-Time PCR Protocol

Follow the Instructions for Use for the PANDAA Ebola kit (Aldatu Biosciences, catalog # 4011096). Refer to Appendix B of the Instructions for Use for detection channel settings for compatible real-time PCR instruments.

Use the real-time PCR template provided by Aldatu Biosciences, available to download from www.aldatu.bio/qpcr-templates.

Data Analysis

Real-Time PCR Analysis Parameters

Refer to Appendix B of the PANDAA Ebola Instructions for Use for instrument-specific analysis settings, including threshold settings for determination of Ct values for the Ebolavirus, Internal Control, and Sample Adequacy Control targets. For technical assistance, contact Aldatu Biosciences Technical Support at support@aldatubio.com.

Expected Results

Label	Virus	Isolate	Expected Result
A	Bundibugyo virus	Butalya-811250	Detected
B	Reston virus	Pennsylvania	Detected
C	Sudan virus	Boniface	Detected
D	Tai Forest virus	Pauléoula-CI	Detected
E	Ebola virus	Makona	Detected
F	Sudan virus	Gulu	Detected
G	Marburg virus	Musoke	Not Detected
H	Ebola virus	Mayinga	Detected
I	Ebola virus	Mayinga	Detected
J	Ebola virus	Mayinga	Detected
Positive Control	Ebola virus	Mayinga	Detected
Negative Control	-	-	Not Detected

Customer and Technical Support

Email: support@aldatubio.com

Address: Aldatu Biosciences
313 Pleasant Street
Watertown, MA 02472
USA

Disclaimers

This Qualification Panel is for Research Use Only and is not intended for use in diagnostic procedures.

Aldatu Biosciences will not be liable for any direct, indirect, consequential, or incidental damage related to or arising from the use, the results of use of this product or document. This User Guide Document is provided “as is” and is subject to being changed, without notice, in future editions. Aldatu Biosciences does not guarantee in any way that you will obtain satisfactory results from using this product as described herein. The only warranty provided is limited product warranty, that is, if this product does not meet the standard performance for control samples, the product will be replaced at no charge.

Trademarks

PANDAA® is a registered trademark of Aldatu Biosciences, Inc. PANDAA EBOV, PANDAA Ebola, and the Aldatu Biosciences logo are trademarks of Aldatu Biosciences, Inc. FAM™ and ROX™ are trademarks of Life Technologies Corporation. ABI Prism® (Applied Biosystems); CFX96™ (Bio-Rad); Cy® (GE Healthcare); LightCycler® (Roche); Rotor- Gene®. All other trademarks that appear in this document are the property of their respective owners.

Patents

The PANDAA technology is covered by US Patent No. 10,100,349 and European Patent Application No. 3052656 owned by the President and Fellows of Harvard College and exclusively licensed to Aldatu Biosciences, Inc.

Copyright

© 2026 Aldatu Biosciences, Inc.

Explanation of Symbols



For Research Use Only



Catalog number



Batch number



Storage temperature range



Use-by date



Contents sufficient for n reactions



Consult Instructions for Use



Contents of the PANDAA® EBOV Qualification Panel



Manufacturer

The logo for ALDATU BIOSCIENCES is displayed in white on a dark blue background. The word "ALDATU" is written in a large, bold, sans-serif font. A horizontal line runs through the middle of the letters "A", "L", "D", "A", and "T". To the left of the "A", there are three curved lines that sweep upwards and to the right, suggesting motion or a biological process. Below "ALDATU", the word "BIOSCIENCES" is written in a smaller, all-caps, sans-serif font with wide letter spacing.

ALDATU
BIOSCIENCES

al • da • tu [ɒl - də - tu]: to become something different